




California Vaccines for Children Program

Digital Data Logger (DDL) Training Resources for CA VFC Providers Receiving a PPHF Grant-funded Data Logger

Logger	Training Resources	Documentation	Software	Technical Support
<p>DeltaTrak FlashLink Data Logger with Glycol Bottle- Model #20936</p> 	<ul style="list-style-type: none"> • Product overview (3 min) 	<ul style="list-style-type: none"> • FlashLink Quick Start Guide • Accuracy does not meet CA VFC requirements if used in the freezer 	<ul style="list-style-type: none"> • Software installation – getting started (2 min) 	<ul style="list-style-type: none"> • Technical support is available by email at techsupport@deltatrak.com or by calling (800) 390-0804. • Support hours are Monday through Friday 7:30 a.m. to 4:30 p.m. Pacific Time.
<p>LogTagVAC-TRED30-KIT Vaccine Monitoring and Storage System</p> 	<ul style="list-style-type: none"> • MircoDAQ training resource for California VFC Providers • LogTag Training California • California DPH Logger Configuration • Instructions to Print Summary Report 	<ul style="list-style-type: none"> • LogTag TRID30 and TRED30 Quick Start Guide • LogTag Analyzer Software Manual • Resolving USB Installation Issues • Meets CA VFC requirements for refrigerator and freezer 	<ul style="list-style-type: none"> • LogTag Analyzer Software (Version 2.5.19) 	<ul style="list-style-type: none"> • Technical support is available by contacting MicroDAQ.com via their Website or by calling (603) 746-5524 and pressing option "1". • Support Hours: Monday through Friday: 9 a.m. to 4:30 p.m. Eastern Time.
<p>Traceable® Excursion-Trac™ Data logging Thermometer #6430</p> 	<ul style="list-style-type: none"> • Traceable-Excursion-Trac Training Videos 	<ul style="list-style-type: none"> • Traceable Excursion-Trac 6430 Instructions • Meets CA VFC requirements for refrigerator and freezer 	<ul style="list-style-type: none"> • The Excursion-Trac does not require software. • A non-encrypted flash drive is needed to download temperature recordings. • Logger operates with an AC adapter (included). A set of 2 AAA batteries function as a backup power source in case of a power failure. 	<ul style="list-style-type: none"> • Technical support is available by e-mail at Sales@Control3.com, or by calling (281) 482-1714. • Support hours: Monday through Friday 8:30 a.m. to 5:00 p.m. Central Time.