



Vaccine Tips

January 2021

If your data logger download looks like this:

	A	B	C	D	E	F
1	ID:	D16634106763202A				
2	Temperature Unit of Measure is Celsius					
3	Date format is MM/DD/YYYY					
4	Alarm Event Data					
5	Probe	Temp	Time Out	Date Out	Time In	Date In
6	1	-13.33	15:08	4/24/2018	15:08	4/24/2018
7	Temperature Data					
8	P1	P2	Time	Date		
9	4	0	11:40	1/7/2019		
10	4.01	0	11:05	1/7/2019		
11	3.94	0	10:30	1/7/2019		
12	3.95	0	9:55	1/7/2019		
13	4.04	0	9:20	1/7/2019		

Then it's time to replace it.

- ✓ Do not recalibrate data loggers that download into an Excel or CSV file.
- ✓ Replace with a data logger that downloads a Summary Report.

Summary Reports will reduce time needed to report temperature excursions by accurately displaying at a glance:

- recorded temperature data since the device was last reset,
- the total time temperatures have been out of range,
- minimum and maximum temperatures, and
- alarm settings.

A properly configured data logger that meets VFC requirements will help ensure your vaccines are stored appropriately. For more information, refer to the California VFC Program [Provider Operations Manual](#).



California Department of Public Health
Immunization Branch
Vaccines for Children (VFC) Program
850 Marina Bay Parkway
Richmond, CA 94804

Toll Free Phone:
877-2GET-VFC (877-243-8832)

Toll Free Fax:
877-FAXX-VFC (877-329-9832)